

2023 PRODUCT CATALOG

TOTAL
● ● ● ● ●
Solution

INHEALTH
TECHNOLOGIES®

CONTENTS

Pre-Operative Education	2
"What to Expect" Brochure	2
Educational Videos	2
Surgical Solutions	3
Tracheoesophageal Puncture Set (Sterile)	3
Voice Prosthesis Placement Surgical Kits (Sterile)	3
Indwelling Voice Prostheses (Sterile)	3
Adjustable Bi-Flanged Fistula Prostheses (Sterile)	3
Post-Operative Care	4
Laryngectomy Tubes	4
Marpac Tracheostomy Tube Holder/Collar	4
Kapitex Trachi-Swab Hi-Vis	4
Electrolarynx	5
Post-Operative Care Kits (For Total Laryngectomy)	6
Next Steps	8
Voice Prostheses (Clinician Placed)	8
Special Order Voice Prostheses (Clinician Placed)	10
Voice Prostheses (Patient Changeable)	12
Heat & Moisture Exchange (HME)	14
Voice Prosthesis Accessories	17
Adhesive Housings	18
Tracheostoma Buttons	19
Tracheostoma Attachments	19
And Beyond	21
Hands-Free Speech Valves	21
Skin Care	22
Stoma Protection and Accessories	23
Rx Only	25
Voice Prosthesis Puncture Kits and Management	26
RENÚ® Voice and RENÚ® Gel	29
Oto-Rhino Foreign Body Remover	30
Nasal Septal Perforation Prostheses	32
Ordering Information	33
New Patient Enrollment Instructions	33
Free Educational Resources	33

InHealth Technologies® offers a TOTAL Solution for your total laryngectomy from pre- to post-op and beyond

TOTAL Solution serves patients at all phases of voice restoration, from products that deliver solutions for pre- and post-operative care and beyond. This continuum of care is enhanced by InHealth Technologies' customer-focused team, including leading clinical specialists, experienced patient consultants and dedicated customer service members.



TOTAL Solution

Pre-Operative Education

- "What to Expect" Brochure
- Educational Videos

PHASE

1

PHASE

2

Surgical Solutions

- Voice Prosthesis Placement Surgical Kits
- Sterile Option Voice Prostheses
- Fistula Prostheses

PHASE

3

Post-Operative Care

- Post-Operative Care Kits
- Laryngectomy Tubes
- Electrolarynx
- "Life After Laryngectomy" Brochure

PHASE

4

Next Steps

- Next Steps Kit
- Voice Prostheses
- HMEs
- Adhesive Housings
- Educational Programs

PHASE

5

And Beyond


- Hands-Free Speech Valve
- Skin Care

● PRE-OPERATIVE EDUCATION

“What to Expect” Brochure

Available online at:

<https://inhealth.com/educational-materials>



Educational Videos

Available at InHealth Technologies®

<https://inhealth.com/videos>



● SURGICAL SOLUTIONS

Blom-Singer® Tracheoesophageal Puncture Set (Sterile)

Designed for secondary placement of a tracheoesophageal prosthesis. Includes a disposable 16 gauge curved introducer needle with a 14 gauge Teflon® sheath, a 26" (66cm) silicone catheter with a tapered reinforced tip attached to a 24" (61cm) Teflon-coated stainless steel leader cable, and a fenestrated surgical drape with adhesive around the aperture.

Ref	Description	Rx
BE 6055	Tracheoesophageal Puncture Set (Sterile)	■



Blom-Singer® Voice Prosthesis Placement Surgical Kits (Sterile)

Voice prosthesis placement for primary or secondary puncture. The Blom-Singer Voice Prosthesis Placement Surgical Kits are designed to be used with Blom-Singer Classic Indwelling Voice Prosthesis, Sterile Series. Sold separately.

Ref	Description	Rx
TP 1000	16 Fr Voice Prosthesis Placement Kit (Sterile)	■
TP 1001	20 Fr Voice Prosthesis Placement Kit (Sterile)	■



Blom-Singer® Classic™ Indwelling Voice Protheses (Sterile)

Designed for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
8mm	IN 1608-SO	IN 2008-SO	■	12mm	IN 1612-SO	IN 2012-SO	■
10mm	IN 1610-SO	IN 2010-SO	■	14mm	IN 1614-SO	IN 2014-SO	■



Blom-Singer® Adjustable Bi-Flanged Fistula Protheses (Sterile)

Designed by Eric D. Blom, Ph.D to control leakage of saliva, food, and drink around a hypopharyngeal fistula. The prosthesis is adjustable in length and comes in three flange sizes, 25mm, 38mm, and 50mm, ensuring a proper fit.

Ref	Description	Rx
FM 1003	25mm (Small) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea	■
FM 1004	38mm (Medium) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea	■
FM 1005	50mm (Large) Adjustable Bi-Flanged Fistula Prosthesis (Available till end of Q4 2023)	■



● POST-OPERATIVE CARE

Blom-Singer® StomaSoft® Laryngectomy Tubes

Blom-Singer StomaSoft Laryngectomy Tubes are flexible, soft, and compatible with HME and other accessories from InHealth Technologies®. They feature reinforced strap holders and are available with or without fenestrations.

Size/Length	I.D./O.D.	Non-Fenestrated	Fenestrated	Rx
8/36	9.5mm/12.0mm	BE 6698	BE 6698F	■
8/55	9.5mm/12.0mm	BE 6699	BE 6699F	■
9/36	10.5mm/13.5mm	BE 6700	BE 6700F	■
9/55	10.5mm/13.5mm	BE 6701	BE 6701F	■
10/36	12.0mm/15.0mm	BE 6702	BE 6702F	■
10/55	12.0mm/15.0mm	BE 6703	BE 6703F	■
12/36	13.5mm/17.0mm	BE 6704	BE 6704F	■
12/55	13.5mm/17.0mm	BE 6705	BE 6705F	■

Blom-Singer® Laryngectomy Tubes (Sterile)

Blom-Singer Laryngectomy Tubes (Sterile) are durable, lightweight silicone tubes that are used to keep the tracheostoma open after laryngectomy. The tubes do not hold an HME.

Size/Length	I.D./O.D.	Non-Fenestrated (Sterile)	Rx
8/36	9.5mm/12.0mm	BE 6398	■
8/55	9.5mm/12.0mm	BE 6399	■
9/36	10.5mm/13.5mm	BE 6400	■
9/55	10.5mm/13.5mm	BE 6401	■
10/36	12.0mm/15.0mm	BE 6402	■
10/55	12.0mm/15.0mm	BE 6403	■
12/36	13.5mm/17.0mm	BE 6404	■
12/55	13.5mm/17.0mm	BE 6405	■

Marpac Tracheostomy Tube Holder/Collar™

Marpac Tracheostomy Tube Holder/Collar is a single use, disposable Comfortable Collar™ used to secure the tracheostomy tube. Latex free. Adult—one size fits all.

Ref	Description
204	Tracheostomy Tube Holder/Collar

Kapitex® Trachi-Swab® Hi-Vis

The Kapitex Trachi-Swab Hi-Vis is indicated for use in patients who have a tracheostoma and use a device to maintain stoma patency. It is intended for cleaning tubes and inner cannulas after removal from the patient. Available in 50 packs.

Ref	Description
TRSWB3050	Trachi-Swab Hi-Vis



Blom-Singer® Electrolarynx EL 1000

An electrolarynx is a speaking aid that produces sound vibrations that can be conducted through the throat and mouth to speak after a laryngectomy.

InHealth Technologies® and Eric D. Blom, Ph.D are pleased to offer a high quality, cost effective electrolarynx for laryngectomees.

Ref	Description
EL 1000	Electrolarynx Digital Speech Aid (One 9-volt batteries included)
EL 1002	Oral Adapter

Oral Adaptor EL 1002



Electrolarynx EL 1000



Blom-Singer® Post-Operative Care Kits (For Total Laryngectomy)

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Post-Operative Care Kits were designed with valued input from clinicians and patients. These kits contain supplies to assist patients in their transition from the hospital to home after a total laryngectomy. They are available in several configurations to meet select needs. All configurations include your choice of StomaSoft Laryngectomy Tube and HMEs as listed below.

The Post-Operative Care Kit is housed in a hard, clear storage case, which also includes information on how to obtain our Next Steps follow-up kit. The Next Steps kit includes additional sample product and educational items.

All kit configurations contain the following patient care supplies:

- 1 each 204 Marpac Tracheostomy Tube Holder/Collar™
- 1 each BE 6048 Blom-Singer® Shower Guard
- 1 each Drawstring Bag
- 1 each Mirror with Stand
- 1 each Halcyon LED Flashlight

REF	Description
PK 2100	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2101	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2102	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2103	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2106	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2107	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2104	
EL 1000	Electrolarynx, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30



PK 2101



● NEXT STEPS

VOICE PROSTHESES (CLINICIAN PLACED)

Blom-Singer® Classic™ Indwelling Voice Prosthesis (Sterile)

Includes sterile prostheses for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
8mm	IN 1608-SO	IN 2008-SO	■	12mm	IN 1612-SO	IN 2012-SO	■
10mm	IN 1610-SO	IN 2010-SO	■	14mm	IN 1614-SO	IN 2014-SO	■



Classic (Sterile)

Blom-Singer® Classic™ Indwelling Voice Prosthesis

Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
4mm	IN 1604-NS	IN 2004-NS	■	14mm	IN 1614-NS	IN 2014-NS	■
6mm	IN 1606-NS	IN 2006-NS	■	16mm	IN 1616-NS	IN 2016-NS	■
8mm	IN 1608-NS	IN 2008-NS	■	18mm	IN 1618-NS	IN 2018-NS	■
10mm	IN 1610-NS	IN 2010-NS	■	20mm	IN 1620-NS	IN 2020-NS	■
12mm	IN 1612-NS	IN 2012-NS	■				




Classic

Blom-Singer® Dual Valve™ Indwelling Voice Protheses

The Blom-Singer Dual Valve Indwelling Voice Prosthesis has two valves embedded with silver oxide as a material preservative. When the primary valve fails, the secondary speech valve is intended to remain viable to help prevent leakage of fluids. The Blom-Singer® Dual Valve Voice Prosthesis is available with standard or large tracheal/esophageal flanges. Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/ cleaning accessories. Clinician placed.

Length	20 Fr	20 Fr Large Flange	Rx
6mm	DV 2006-NS	DV 2006-LF	■
8mm	DV 2008-NS	DV 2008-LF	■
10mm	DV 2010-NS	DV 2010-LF	■

Length	20 Fr	20 Fr Large Flange	Rx
12mm	DV 2012-NS	DV 2012-LF	■
14mm	DV 2014-NS	DV 2014-LF	■



Blom-Singer® Advantage® Indwelling Voice Protheses

The Blom-Singer Advantage Indwelling Voice Prosthesis is embedded with silver oxide as a material preservative.* Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories.

HARD VALVE ASSEMBLY

Length	20 Fr	Rx
4mm	AD 2004-NS	■
6mm	AD 2006-NS	■
8mm	AD 2008-NS	■
10mm	AD 2010-NS	■
12mm	AD 2012-NS	■
14mm	AD 2014-NS	■

SOFT VALVE ASSEMBLY

Length	16 Fr	20 Fr	Rx
4mm	AS 1604-NS	AS 2004-NS	■
6mm	AS 1606-NS	AS 2006-NS	■
8mm	AS 1608-NS	AS 2008-NS	■
10mm	AS 1610-NS	AS 2010-NS	■
12mm	AS 1612-NS	AS 2012-NS	■
14mm	AS 1614-NS	AS 2014-NS	■



*Silver oxide, a known anti-microbial substance, has been incorporated into the silicone flap valve of the Blom-Singer Advantage and Dual Valve Indwelling Voice Prothese as a material preservative.

SPECIAL ORDER VOICE PROSTHESES (CLINICIAN PLACED)


—DEVELOPED BY ERIC D. BLOM, PH.D

Unique problems call for unique solutions. InHealth Technologies® offers several voice prostheses to address those needs.

Special Length

Available in one millimeter incremental lengths for a more precise TE puncture fit.


	5mm	7mm	9mm	Rx
16 Fr	IN 1605-SL	IN 1607-SL	IN 1609-SL	■
20 Fr	IN 2005-SL	IN 2007-SL	IN 2009-SL	■



Increased Resistance

Designed for patients who experience intermittent valve opening with fluctuations in esophageal pressures, resulting in air ingestion and/or leakage.


	6mm	8mm	10mm	12mm	Rx
16 Fr	IN 1606-IR	IN 1608-IR	IN 1610-IR	IN 1612-IR	■
20 Fr	IN 2006-IR	IN 2008-IR	IN 2010-IR	IN 2012-IR	■



Large Esophageal Flange


Provides a seal to help manage leakage around a voice prosthesis. Also, it enables downsizing from a wider diameter prosthesis to a smaller diameter. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LEF	IN 1606-LEF	IN 1608-LEF	IN 1610-LEF	IN 1612-LEF	IN 1614-LEF	■
20 Fr	IN 2004-LEF	IN 2006-LEF	IN 2008-LEF	IN 2010-LEF	IN 2012-LEF	IN 2014-LEF	■



LARGE ESOPHAGEAL FLANGE (SPECIAL LENGTH)

	5mm	7mm	Rx
16 Fr	IN 1605-LESL	IN 1607-LESL	■
20 Fr	IN 2005-LESL	IN 2007-LESL	■



Large Esophageal and Tracheal Flanges

Provides a double seal in patients with enlarged puncture tracts and tear-drop or irregularly shaped punctures. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LF	IN 1606-LF	IN 1608-LF	IN 1610-LF	IN 1612-LF	IN 1614-LF	■
20 Fr	IN 2004-LF	IN 2006-LF	IN 2008-LF	IN 2010-LF	IN 2012-LF	IN 2014-LF	■




LARGE ESOPHAGEAL FLANGE (INCREASED RESISTANCE)

	4mm	6mm	8mm	Rx
16 Fr	IN 1604-LEIR	IN 1606-LEIR	IN 1608-LEIR	■
20 Fr	IN 2004-LEIR	IN 2006-LEIR	IN 2008-LEIR	■

TEP Occluders

The Blom-Singer® TEP (Tracheoesophageal Puncture) Occluder is a non-valved silicone device designed for placement in, and maintenance of, the TEP following total laryngectomy. Not a voice prosthesis.

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LETO	IN 1606-LETO	IN 1608-LETO	IN 1610-LETO	IN 1612-LETO	IN 1614-LETO	■
20 Fr	IN 2004-LETO	IN 2006-LETO	IN 2008-LETO	IN 2010-LETO	IN 2012-LETO	IN 2014-LETO	■



VOICE PROSTHESES (PATIENT CHANGEABLE)

Blom-Singer® Duckbill Voice Prostheses

The Blom-Singer® Duckbill is a patient changeable, slit valve prosthesis.

Length	16 Fr	Rx	Length	16 Fr	Rx
6mm	DB 16-006	■	12mm	DB 16-012	■
8mm	DB 16-008	■	14mm	BE 6001	■
10mm	BE 6000	■	18mm	BE 6002	■



Blom-Singer® Low Pressure Voice Prostheses

The Blom-Singer® Low Pressure is a patient changeable, one-way valve prosthesis.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
6mm	BE 6009	BE 2014	■	18mm	BE 6012	BE 2026	■
8mm	LP 16-008	LP 20-008	■	22mm	BE 6013	BE 2030	■
10mm	BE 6010	BE 2018	■	25mm	BE 6014	BE 2033	■
12mm	LP 16-012	LP 20-012	■	28mm	BE 6015	—	■
14mm	BE 6011	BE 2022	■				



Blom-Singer® Special Order Low Pressure Voice Prostheses

The Blom-Singer® Special Order Low Pressure is a patient changeable, one-way valve prosthesis.

Length	16 Fr Special Increased Resistance	16 Fr Special Length	20 Fr Special Length	20 Fr Special Increased Resistance	Rx
4mm	—	SPECIAL 11-902-06	SPECIAL 11-251-06	—	■
6mm	LP 16-006-IR	—	—	LP 20-006-IR	■
8mm	LP 16-008-IR	—	—	LP 20-008-IR	■
10mm	LP 16-010-IR	—	—	LP 20-010-IR	■
12mm	LP 16-012-IR	—	—	LP 20-012-IR	■
14mm	LP 16-014-IR	—	—	LP 20-014-IR	■
16mm	—	SPECIAL 11-902-04	SPECIAL 11-251-04	—	■
18mm	LP 16-018-IR	—	—	—	■
20mm	—	SPECIAL 11-902-05	SPECIAL 11-251-05	—	■
22mm	LP 16-022-IR	—	—	—	■
25mm	LP 16-025-IR	—	—	—	■
28mm	LP 16-028-IR	—	—	—	■

Blom-Singer® Voice Prosthesis Sizers

The Blom-Singer Voice Prosthesis Sizer is a device used by clinicians to measure the length of the tracheoesophageal puncture to assist in selection of the appropriate voice prosthesis.

Ref	Description	Rx
BE 6036	Voice Prosthesis Sizer with 3 inserters (for 16 Fr voice prostheses)	■
BE 2036	Voice Prosthesis Sizer with 3 inserters (for 20 Fr voice prostheses)	■
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea	■

Blom-Singer® Dilator-Sizers

The Blom-Singer Dilator-Sizer combines two tasks into one innovative product: sizing and dilation. In one step, the clinician can dilate and measure the tracheoesophageal puncture during voice prosthesis placement procedures.

Ref	Description	Rx
BE 2051	18 Fr Dilator-Sizer (for 16 Fr prostheses)	■
BE 2052	22 Fr Dilator-Sizer (for 20 Fr prostheses)	■

VOICE PROSTHESIS ACCESSORIES

Blom-Singer® Replacement Gel Caps

Use of the Blom-Singer Replacement Gel Caps provides a smooth, rounded shape to the esophageal flange end of the voice prosthesis, enabling easy insertion into the tracheoesophageal puncture.

Ref	Description	Rx
BE 3190	Replacement Gel Caps, 16 Fr, 90 ea	■
BE 3195*	Replacement Gel Caps, 18 Fr, 90 ea	■
BE 3290	Replacement Gel Caps, 20 Fr, 90 ea	■
BE 3295*	Replacement Gel Caps, 22 Fr, 60 ea	■

*BE 3195 and BE 3295 are not compatible with LP 3100.

Blom-Singer® Flushing Device and Cleaning Brushes

For routine cleaning of a voice prosthesis as recommended by your clinician.

Ref	Description
IN 4050	Flushing Device, pack of 3
IN 4064*	Cleaning Brush, 16/20 Fr, Small, 4mm–8mm, 3 ea
IN 4065*	Cleaning Brush, 16/20 Fr, Medium, 9mm–16mm, 3 ea
IN 4066*	Cleaning Brush, 16/20 Fr, Large, 18mm–22mm, 3 ea

*Fits indwelling style voice prostheses only.

Blom-Singer® Plug Inserts

The Blom-Singer Plug Insert is designed for temporary management of leakage through a voice prosthesis. Not compatible with patient changeable voice prostheses.

Ref	Description	Rx
IN 4069	16 Fr Plug Insert (Advantage and Classic, 6mm to 20mm)	■
IN 4070	20 Fr Plug Insert (Advantage and Classic, 6mm to 20mm)	■



Blom-Singer® HME All Day/All Night Advantages around the clock



ClassicFlow® HME



EasyFlow® HME



Day&Night™
ClassicFlow® HME



Day&Night™
EasyFlow® HME

- 24-hour use
- EasyTouch™ speech button—push to speak
- Compatible with a wide variety of accessories

HEAT & MOISTURE EXCHANGE (HME)

Research shows that long-term HME use results in less mucus¹, less coughing², and decreased risk of shortness of breath, fatigue, sleep problems, anxiety, and depression^{1,2}.

Blom-Singer® SpeakFree® Hands-Free Valve with ClassicFlow® HME

Blom-Singer SpeakFree HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. ClassicFlow is designed for increased moisture retention and resistance. The ClassicFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1090EZ	SpeakFree Hands-Free Valve with ClassicFlow HME, 30 per pack

Blom-Singer® SpeakFree® Hands-Free Valve with EasyFlow® HME

Blom-Singer SpeakFree HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. With EasyFlow, the foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs. The EasyFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1090EF	SpeakFree Hands-Free Valve with EasyFlow HME, 30 per pack

Blom-Singer® Day&Night™ ClassicFlow® HME

The Blom-Singer Day&Night ClassicFlow can be worn all day and night including bedtime. The ClassicFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1091	Day&Night ClassicFlow HME Cartridge, pack of 30 (24-hour use)

Blom-Singer® Day&Night™ EasyFlow® HME

The Blom-Singer Day&Night EasyFlow allows you to breathe more freely to match your activity level and pulmonary needs. The EasyFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1092	Day&Night EasyFlow HME Cartridge, pack of 30 (24-hour use)



NEW

ClassicFlow



NEW

EasyFlow



NEW

ClassicFlow



NEW

EasyFlow

1. Ackerstaff AH, Hilgers FJ, et al. Improvements in respiratory and psychosocial functioning following total laryngectomy by the use of a heat and moisture exchanger. *Annals of Otolaryngology, Rhinology & Laryngology*. 1993;102(11):878–883.
2. Ackerstaff AH, Hilgers FJ, et al. Heat and moisture exchangers as a treatment option in the post-operative rehabilitation of laryngectomized patients. *Clinical Otolaryngology*. 1995;20:504–509.

Blom-Singer® ClassicFlow® HME (Available till end of 2023)

The Blom-Singer ClassicFlow HME is designed for increased moisture retention and airflow resistance.

Ref	Description
BE 1055EZ	ClassicFlow HME Cartridge, pack of 30 (24-hour use)



Blom-Singer® EasyFlow® HME (Available till end of 2023)

The foam filter of Blom-Singer EasyFlow allows you to breathe more freely to match your activity level and pulmonary needs. The EasyFlow has a non-reflective beige tint.

Ref	Description
BE 1055EF	EasyFlow HME Cartridge, pack of 30 (24-hour use)



Blom-Singer® HumidiFilter® Holder (Available till end of 2023)

The Blom-Singer HumidiFilter Holder is a reusable, lightweight plastic holder that holds a HumidiFilter foam filter (assembly required).

Ref	Description
BE 1060	HumidiFilter Holder
BE 1080	HumidiFilter Foam Filters, pack of 30 (Not for 24-hour use)





ADHESIVE HOUSINGS

InHealth Technologies® presents the FitSeries—adhesive housings designed to meet your skin and stoma needs after a total laryngectomy. Available in two styles, three shapes each. Compatible with HME cartridges, SpeakFree HME Hands-Free, ATSV and ATSV II.

Blom-Singer® AccuFit® Adhesive Housings

The Blom-Singer AccuFit Adhesive Housing is designed to be our strongest adhesive option while still maintaining flexibility. It has increased elasticity and is less prone to tearing. Available in three shapes: round, oval, and oval extra.

Ref	Description
BE 6082	AccuFit Adhesive Housing, Round, pack of 30
BE 6083	AccuFit Adhesive Housing, Oval, pack of 30
BE 6084	AccuFit Adhesive Housing, Oval Extra, pack of 30

Blom-Singer® HydroFit® Adhesive Housings


The Blom-Singer HydroFit Adhesive Housing is a hydrocolloid adhesive baseplate made with a protective gel-like layer and a moisture-absorbent adhesive. The HydroFit is gentle on the skin and flexible, but maintains strong adhesion. Available in three shapes: round, oval, and oval extra.

Ref	Description
BE 6085	Hydrofit Adhesive Housing, Round, pack of 30
BE 6086	Hydrofit Adhesive Housing, Oval, pack of 30
BE 6087	Hydrofit Adhesive Housing, Oval Extra, pack of 30

Blom-Singer® TruSeal® Contour™ Low Profile Adhesive Housings (Available till end of 2023)

The Blom-Singer TruSeal Contour Low Profile Adhesive Housing is a lightweight, flexible material made to fit various peristomal anatomies. It has a lower profile from the standard housing, reducing its height and overall profile. Available in three shapes: round, oval, and standard.

Ref	Description
BE 6074	TruSeal Contour Low Profile Adhesive Housing, Round, pack of 30
BE 6075	TruSeal Contour Low Profile Adhesive Housing, Oval, pack of 30
BE 6076	TruSeal Contour Low Profile Adhesive Housing, Standard, pack of 30



TRACHEOSTOMA BUTTONS

Barton-Mayo™ Tracheostoma Buttons*

Developed by Medical Innovations International, Inc. in association with the Mayo Clinic to provide tracheostoma valve retention for near-total laryngectomy and total laryngectomy patients.

Ref	Size	O.D.	Length
BMB091	9	13.5mm	15mm, short
BMB092	9	13.5mm	22mm, regular
BMB093	9	13.5mm	30mm, long
BMB101	10	15mm	15mm, short
BMB102	10	15mm	22mm, regular
BMB103	10	15mm	30mm, long
BMB121	12	17mm	15mm, short
BMB122	12	17mm	22mm, regular
BMB123	12	17mm	30mm, long
BMB141	14	19mm	15mm, short
BMB142	14	19mm	22mm, regular
BMB143	14	19mm	30mm, long
BMBKIT Sizing Kit (all 12 sizes)			
BMBKIT-SHORT Patient Sizing Kit-Short (4 short sizes only)			

*Only available for sale in the US.



TRACHEOSTOMA ATTACHMENTS

Blom-Singer® Tracheostoma Valve Housings

The Tracheostoma Valve Housing is a reusable, flexible housing that attaches the ATSV II, HME cartridge or Shower Guard over the stoma. It is applied with double sided tape discs. It is available in two sizes, standard and large.

Ref	Description
BE 6038	Tracheostoma Valve Housing, standard, 1 ea
BE 6039	Tracheostoma Valve Housing, large, 1 ea



Blom-Singer® Adhesive Discs

A variety of double-sided tape discs are available to provide attachment of the tracheostoma valve housing.

Ref	Description	Ref	Description
BE 6041	Regular Tape Discs, standard, pack of 30	BE 6044	Thin Tape Discs, large, pack of 30
BE 6042	Regular Tape Discs, large, pack of 30	BE 6034	Heavy Duty Tape Discs, standard, pack of 30
BE 6043	Thin Tape Discs, standard, pack of 30	BE 6035	Heavy Duty Tape Discs, large, pack of 30





● AND BEYOND

HANDS-FREE SPEECH VALVE

Blom-Singer® SpeakFree® Hands-Free Valve with ClassicFlow® HME

Blom-Singer SpeakFree HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. ClassicFlow is designed for increased moisture retention and resistance. The ClassicFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1090EZ	SpeakFree Hands-Free Valve with ClassicFlow HME, 30 per pack



ClassicFlow

NEW

Blom-Singer® SpeakFree® Hands-Free Valve with EasyFlow® HME

Blom-Singer SpeakFree HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. With EasyFlow, the foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs. The EasyFlow is translucent and non-reflective for a natural look.

Ref	Description
BE 1090EF	SpeakFree Hands-Free Valve with EasyFlow HME, 30 per pack



EasyFlow

NEW

Blom-Singer® Adjustable Tracheostoma Valve II (ATSV II)

The Adjustable Tracheostoma Valve II (ATSV II) consists of two pieces: a Diaphragm/Faceplate and Body. The assembled two-piece valve inserts into the tracheostoma adhesive housing that is applied to the skin surrounding the stoma. A replacement Diaphragm/Faceplate is available separately.

Ref	Description	Rx
BE 8025H	ATSV II Valve (Body & Diaphragm/Faceplate), 1 kit	■
BE 8026H	ATSV II Replacement Diaphragm/Faceplate, 1 ea	■



BE 8025H

Blom-Singer® ATSV II HumidiFilter® Cap and 7 Foam Filters

The ATSV II HumidiFilter Cap is a twist-on plastic cap that fits over the ATSV II. Comes with seven HumidiFilter Foam Filters. Additional replacement filters sold separately. (For use with ATSV II only.)

Ref	Description
BE 1010	ATSV II HumidiFilter Cap and 7 Foam Filters, 1 kit
BE 1030	ATSV II HumidiFilter Foam Filters, pack of 30



SKIN CARE

Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes

Skin-Prep is a liquid film-forming dressing that, upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films.

Ref	Description
59420425	Skin-Prep Liquid Film-Forming Protective Barrier Wipes, box of 50



Skin Tac™ Adhesive Barrier Wipes

Skin Tac liquid adhesive barrier wipes are latex-free, hypo-allergenic, clear, non-rubber adhesive barrier wipe.

Ref	Description
MS407W*	Adhesive Barrier Wipes, box of 50

*May not be available in all international markets.



NFS Adhesive (Silicone Adhesive)

PSA Pressure Sensitive Adhesive, silicone based adhesive for use with all silicone prosthetic devices.

Ref	Description
B-530	NFS Adhesive, 1 oz



Remove™ Adhesive Remover

Gently cleans hydrocolloid, acrylic and rubber-based adhesive residues from the skin without irritation, discomfort or trauma.

Ref	Description
59403125	Adhesive Remover, box of 50



STOMA PROTECTION AND ACCESSORIES

ADDvox® Stoma Cover Mock Turtlenecks (Only available in the US)

The ADDvox Stoma Cover Mock Turtleneck has inner layers of material which helps keep debris from entering the stoma. They slip under a shirt like a fashion dickey and have an outer layer of soft 100% cotton knit. Each cover is 9-1/2" long x 8-3/4" wide. Velcro fastener adjusts to fit necks of 12-1/2" to 22". Each color sold separately.

Ref	Description	Color	Quantity
BE 6227-R2	Stoma Cover Mock Turtleneck	White	1 ea
BE 6228-R2	Stoma Cover Mock Turtleneck	Black	1 ea
BE 6229-R2	Stoma Cover Mock Turtleneck	Navy	1 ea



ADDvox® Stoma Scarves (Only available in the US)

The ADDvox Stoma Cover Scarf has inside layers to help filter out dust, dirt and pollen, is completely washable, and needs little care. This attractive neckline scarf looks just like a regular scarf or ascot and is available in a wide range of colors for men or women. All scarves have Velcro closures and adjust to fit necks from 17-3/4" to 22" (can be made smaller with tucks or stitching). Each color sold separately.

Ref	Description	Color	Quantity
BE 6232-R2	Stoma Cover Scarf	Black	1 ea
BE 6233-R2	Stoma Cover Scarf	White	1 ea
BE 6234-R2	Stoma Cover Scarf	Navy	1 ea



ADDvox® Cloth Bib-Style Stoma Covers (Only available in the US)

The ADDvox soft, long-lasting bib protects the tracheostoma and lets you breathe more easily. This "bib-style" Stoma Cover will look crisp and clean after repeated laundering. Five soft layers keep dust and dirt from reaching the trachea. Fits necks up to 22". Velcro closure. Each color sold separately.

Ref	Description	Color	Quantity
BE 6224-R2	Cloth Bib-Style Stoma Cover	White	1 ea
BE 6225-R2	Cloth Bib-Style Stoma Cover	Beige	1 ea
BE 6226-R2	Cloth Bib-Style Stoma Cover	Navy	1 ea



ADDvox® Silicone Shower Shield (Only available in the US)

The one-piece ADDvox Silicone Shower Shield is worn during showering to deflect the water away from your stoma. Extremely soft and pliable, it creates a soft, comfortable and conforming seal along the irregular contours of your neck, and vents on the underside allowing for airflow. Made of medical grade silicone, this latex-free Silicone Shower Shield is comfortable, durable, and not affected by temperature, ozone, sunlight, soap or detergents.



Ref	Description
BE 6079-R2	Silicone Shower Shield, 1 ea

Blom-Singer® Shower Guard

The Blom-Singer Shower Guard is a one-piece, right-angle plastic shower guard that comes with a Tracheostoma Standard Adhesive Housing and 30 double-sided Adhesive Tape Discs. Also compatible with the Blom-Singer StomaSoft® Laryngectomy Tube.



Ref	Description
BE 6048	Shower Guard, 1 ea

Blom-Singer® Foam Stoma Protector

A simple and temporary (disposable) solution to help prevent debris from entering the stoma. Adhesive strip on one edge allows for easy attachment, access to the stoma to occlude for speech with a voice prosthesis, or for coughing. Measurement listed is approximate.



Ref	Description
BE 6220	Foam Stoma Protector, Tan, Large 2.7" L x 2.7" W (7cm x 7cm), pack of 30

ADDvox® Stoma Filters with Microporous Adhesive

ADDvox Stoma Filters offer protection against dust, dirt and pollen. Designed for single use with a hypoallergenic, non-irritating, self-adhesive microporous strip that is comfortable and easy to apply. Not recommended for cannula users. 2-1/4" wide x 2-1/2" high x 1/8" thick. Available in off-white and skin tone.



Ref	Description
BE 6222-R2	Stoma Filters, Skin tone, pack of 60
BE 6223-R2	Stoma Filters, Off-white, pack of 60



 ONLY

VOICE PROSTHESIS PUNCTURE KITS AND MANAGEMENT


Blom-Singer® Tracheoesophageal Puncture Set (Sterile)

Ref	Description	Rx
BE 6055	Tracheoesophageal Puncture Set	■



Blom-Singer® Voice Prosthesis Placement Surgical Kits (Sterile)

Ref	Description	Rx
TP 1000	16 Fr Voice Prosthesis Placement Kit (Sterile)	■
TP 1001	20 Fr Voice Prosthesis Placement Surgical Kit (VPP)	■



Blom-Singer® Voice Prosthesis Sizers

Ref	Description	Rx
BE 6036	Voice Prosthesis Sizer with 3 inserters (for 16 Fr voice prostheses)	■
BE 2036	Voice Prosthesis Sizer with 3 inserters (for 20 Fr voice prostheses)	■
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea	■



Blom-Singer® Dilator-Sizers

Ref	Description	Rx
BE 2051	18 Fr Dilator-Sizer (for 16 Fr prostheses)	■
BE 2052	22 Fr Dilator-Sizer (for 20 Fr prostheses)	■



Blom-Singer® Replacement Gel Caps

Ref	Description
BE 3190	Replacement Gel Caps, 16 Fr, 90 ea
BE 3195*	Replacement Gel Caps, 18 Fr, 90 ea
BE 3290	Replacement Gel Caps, 20 Fr, 90 ea
BE 3295	Replacement Gel Caps, 22 Fr, 60 ea

*BE 3195 is not compatible with LP 3100.



Blom-Singer® Tracheoesophageal Puncture Dilators

Ref	Description	Rx
BE 6050	18 Fr Tracheoesophageal Puncture Dilator (for 16 Fr prostheses)	■
BE 2050	22 Fr Tracheoesophageal Puncture Dilator (for 20 Fr prostheses)	■



Blom-Singer® Post-Operative Care Kits (For Total Laryngectomy)

All kit configurations contain the following patient care supplies:

- 1 each 204 Marpac Tracheostomy Tube Holder/Collar™
- 1 each BE 6048 Blom-Singer® Shower Guard
- 1 each Drawstring Bag
- 1 each Mirror with Stand
- 1 each Halcyon LED Flashlight


Ref	Description
PK 2100	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2101	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2102	
BE 6701	StomaSoft Laryngectomy Tube 9/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2103	
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2106	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30
PK 2107	
BE 6705	StomaSoft Laryngectomy Tube 12/55 Non-fenestrated, 1 ea
BE 1055EF	HME Cartridge EasyFlow Heat and Moisture Exchange, pack of 30
PK 2104	
EL 1000	Electrolarynx, 1 ea
BE 6703	StomaSoft Laryngectomy Tube 10/55 Non-fenestrated, 1 ea
BE 1055EZ	HME Cartridge ClassicFlow Heat and Moisture Exchange, pack of 30



PK 2101

Blom-Singer® Adjustable Bi-Flanged Fistula Protheses (Sterile)

Ref	Description
FM 1003	25mm (Small) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1004	38mm (Medium) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1005	50mm (Large) Adjustable Bi-Flanged Fistula Prosthesis



Manometer with Manometer Adaptor

A tool to teach patients how to self-regulate intratracheal air pressure. Proper airflow control may allow for better voice quality, less fatigue, and longer duration of the tracheostoma adhesive housing seal.

Ref	Description
BE 8000	Manometer Kit includes 1 manometer, 1 manometer adaptor and 1 adhesive housing
BE 8020	Replacement Manometer Adaptor includes 5 manometer adaptors and 5 adhesive housings



Blom-Singer® Insufflation Test Set

The Blom-Singer Insufflation Test Set is designed to aid the clinician in subjectively predicting the ability of the laryngectomized patient to produce tracheoesophageal voice.

Ref	Description
BE 6040	Insufflation Test Set, box of 6



The Choice for Long- or Short-Term Augmentation

RENÚ® VOICE AND RENÚ® GEL

RENÚ Voice is for long-term augmentation:

- A buffered hydrogel of synthetic, proven-biocompatible calcium hydroxylapatite (CaHA) particles 25 to 45 microns in size for consistent and easy injection.¹

RENÚ Gel is for short-term augmentation:

- A buffered hydrogel of sodium carboxymethylcellulose² (NaCMC) and glycerin.
- Resorbs slowly in vivo, providing a scaffold for local tissue infiltration.³

RENÚ® Voice: Calcium Hydroxylapatite Implant

RENÚ Voice is an injectable vocal fold implant intended for long-term augmentation.

Ref	Description
08-015-04-V04	RENÚ Voice Calcium Hydroxylapatite Implant, 1.5cc
08-015-04-V01	RENÚ Voice Calcium Hydroxylapatite Implant, 1.5cc, CE

RENÚ® Gel: Hydrogel Implant

RENÚ Gel is an injectable vocal fold implant intended for short-term augmentation.

Ref	Description
08-015-00-V04	RENÚ Gel Hydrogel Implant, 1.5cc

1. Data on file, InHealth Technologies®.

2. Mallur PS, et al. Safety and efficacy of carboxymethylcellulose in the treatment of glottic insufficiency. *The Laryngoscope*. 2012;122:322–326. <https://doi.org/10.1002/lary.21930>.

3. Hughes RGM, Morrison M. Vocal cord medialization by transcutaneous injection of calcium hydroxylapatite. *Journal of Voice*. 2005;19(4):674–678.

RENÚ® VOICE





RENÚ® GEL



Katz Extractor[®]

OTO-RHINO FOREIGN BODY REMOVER

The easy way out





Katz Extractor® Oto-Rhino Foreign Body Remover

Designed in conjunction with David L. Katz, MD, MPH, the Katz Extractor Oto-Rhino Foreign Body Remover is intended for extraction of impacted foreign bodies in the nasal passage. Single use only. This device may also be used in the external ear canal, based upon clinical judgment. Covered by one or more of the following patents and other patents pending: US#5,454,817. Katz Extractor® is a registered trademark of Freudenberg Medical, LLC.

Ref	Description
FB 5000	Oto-Rhino Foreign Body Remover, box of 5




Nasal Septal Perforation Prostheses


The Blom-Singer® Nasal Septal Perforation Prosthesis is a two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

OVAL


Ref	Description
SP 2112	Nasal Septal Perforation Prosthesis Oval 11 x 13 x 2mm
SP 2113	Nasal Septal Perforation Prosthesis Oval 11 x 13 x 3mm
SP 2132	Nasal Septal Perforation Prosthesis Oval 13 x 17 x 2mm
SP 2133	Nasal Septal Perforation Prosthesis Oval 13 x 17 x 3mm
SP 2152	Nasal Septal Perforation Prosthesis Oval 15 x 21 x 2mm
SP 2153	Nasal Septal Perforation Prosthesis Oval 15 x 21 x 3mm




Oval 11 x 13mm



Oval 13 x 17mm




Oval 15 x 21mm




ROUND


Ref	Description
SP 1052	Nasal Septal Perforation Prosthesis Round 5 x 2mm
SP 1053	Nasal Septal Perforation Prosthesis Round 5 x 3mm
SP 1072	Nasal Septal Perforation Prosthesis Round 7 x 2mm
SP 1073	Nasal Septal Perforation Prosthesis Round 7 x 3mm
SP 1102	Nasal Septal Perforation Prosthesis Round 10 x 2mm
SP 1103	Nasal Septal Perforation Prosthesis Round 10 x 3mm




Round 5mm



Round 7mm



Round 10mm



ORDERING INFORMATION

UNITED STATES

Phone: Toll-Free: 800.477.5969

Fax: Toll-Free: 888.371.1530

E-mail: order@inhealth.com

Prescription E-mail: ihtrx@inhealth.com

Order Online: www.inhealth.com

Mail:

InHealth Technologies®
1110 Mark Ave.
Carpinteria, CA 93013 • USA

CANADA

Phone: Toll-Free: 1.888.437.4123

Fax: Toll-Free: 1.226.641.5285

E-mail: canadaorder@inhealth.com

Website: www.inhealth.com/canada

Mail:

InHealth Technologies®
65 Spruce Street
Tillsonburg ON N4G 5C4

**Rx: Products
with the Rx
notation require
a prescription
from a physician
(US patients).**

NEW PATIENT ENROLLMENT INSTRUCTIONS

Enrolling as a new patient? It's easy! Here is what to do:

1. Obtain a **prescription** from your provider.
2. Fill out the **Patient Enrollment Form**.
3. Send back to InHealth Technologies® with a copy of your insurance card (front and back). You can send via email: order@inhealth.com or fax: **888.371.1530**
4. Place your order! Please remember that a prescription alone will not automatically generate an order.

Need insurance or Medicare help? Contact us!

Phone: **800.477.5969 option 3**

Email: medicare@inhealth.com

While you don't have to fill anything out on the **Patient Bill of Rights** form, please make sure you review it. We want you to know your rights.


All forms are available under the **Resources** tab on the InHealth Technologies® website: www.inhealth.com

FREE EDUCATIONAL RESOURCES

- Available on www.inhealth.com under the **Resources** tab.
- Visit InHealth Technologies® video channel at <https://inhealth.com/videos> to view our videos

In 1978, Drs. Singer and Blom pioneered the methods of tracheoesophageal puncture and valved silicone voice prostheses that have been the international standard for voice restoration for over 40 years, allowing thousands of patients worldwide to regain their ability to speak. As a leading supplier of advanced voice restoration systems, our goal at InHealth Technologies® is to help laryngectomees regain and preserve a vital part of their lives—their voice. We are proud to continue the Blom-Singer® legacy today with clinical support, input, and leadership.

Blom-Singer®





STAY CONNECTED...

Speakers Club™

Patients: Please visit us at <https://inhealth.com>

We created this digital E-News for you! **Speakers Club E-News** will provide educational laryngectomy and voice prosthesis information.

ProNET™

Professionals: Please visit us at <https://inhealth.com>

We created this digital E-News specifically for professionals. **ProNET E-News** will deliver useful laryngectomy and voice prosthesis information for you and your practice, including material you may want to share with your patients.

INHEALTH
TECHNOLOGIES®

United States
1110 Mark Ave., Carpinteria, CA 93013
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

International
Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: international@inhealth.com

Patents: <https://inhealth.com/media/documents/pdf/Patents.pdf>
Terms & Conditions: https://inhealth.com/media/documents/pdf/Terms_and_Conditions.pdf
©2023 InHealth Technologies®—Manufactured by Freudenberg Medical, LLC (222002.01)

Blom-Singer®

The voice of experience since 1978



Order online at www.inhealth.com