

GooseNeck®
Snare

We'll pull you
through
Snare and Microsnare



GooseNeck® Snare and Microsnare

Amplatz GooseNeck® Snare Kit (kit includes Snare, Snare Catheter)

Model No.	Loop Diameter (mm)	Snare Length (cm)	Catheter Size (F)	Catheter Length (cm)
GN500	5	120	4	102
GN1000	10	120	4	102
GN1500	15	120	6	102
GN2000	20	120	6	102
GN2500	25	120	6	102
GN3000	30	120	6	102
GN3500	35	120	6	102

Amplatz GooseNeck® Snare Kit - Petite (kit includes Snare, Snare Catheter)

Model No.	Loop Diameter (mm)	Snare Length (cm)	Catheter Size (F)	Catheter Length (cm)
GN1001	10	65	4	48
GN2501	25	65	6	48

Amplatz GooseNeck® Microsnare Kit (kit includes Snare, Microcatheter)

Model No.	Loop Diameter (mm)	Snare Length (cm)	Catheter Size (F)	Catheter Length (cm)
SK200	2	175	2.3 – 3	150
SK201	2	200	2.3 – 3	175
SK400	4	175	2.3 – 3	150
SK401	4	200	2.3 – 3	175
SK700	7	175	2.3 – 3	150
SK701	7	200	2.3 – 3	175

Snare Replacement Catheter

Model No.	Catheter Size (F)	Catheter Length (cm)
MC4000	4	102
MC4001	4	48
MC6000	6	102
MC6001	6	48

Specifications Nominal

Wide Variety of Sizes for Optimal Fit

- Proper sizing of the snare prior to procedure ensures straight-forward retrieval in target vessel and less manipulation time within the vessel

Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

ev3 Corporate World Headquarters / Peripheral Vascular PH +1 763 398 7000 Cust Srvc +1 800 716 6700
ev3 Neurovascular PH +1 949 837 3700
ev3 Europe International Headquarters PH +33 156 88 59 10


Gold Tungsten Coils for Enhanced Visualization

- Gold tungsten loops provide excellent radiopacity for improved procedural precision



True 90° Snare Loop

- Snare loop forms a true 90° angle
- Device remains coaxial to the lumen for proper insertion and successful retrieval or manipulation of atraumatic foreign bodies



Nitinol Shaft for Durability and Kink Resistance

- Super-elastic and shape memory properties of nitinol provide kink resistance